

Multiplex detection of KRAS mutations in colorectal cancer FFPE and cfDNA samples using droplet digital PCR

BIO-RAD

Dawne N Shelton¹, Helen Huang², Wei Yang¹, Jennifer R Berman¹, Samantha Cooper¹, Eli Hefner¹, Filip Janku², John F Regan¹


¹Digital Biology Center, Bio-Rad Laboratories, 5731 W. Las Positas Blvd, Pleasanton, CA, 94588

²MD Anderson Cancer Center, UT Dept. of Investigational Cancer Therapeutics, Houston, TX, 77030

Introduction

Targeted therapies in many cancers have allowed unprecedented progress in the treatment of disease. However, routine implementation of genomic testing is limited due to: 1) difficulties in detection of mutational loads below 5%, 2) limited amounts of sample (pg-ng range) per biological specimen, 3) diagnostic turnaround time and workflow, and 4) cost. KRAS is mutated in approximately 40% of colorectal cancers, and KRAS mutations are predictive of a negative response to α EGFR therapy. To optimize therapy strategies for personalized care, it is therefore critical to rapidly screen patient samples for the presence of multiple KRAS mutations. We have developed a multiplexing strategy to screen clinically-actionable KRAS mutations using digital PCR. No pre-amplification step is required. This sensitive and inexpensive method reduces the risk of contamination and can be easily implemented in molecular diagnostic laboratories for rapid, routine screening of cancer patients.

Materials and Methods



- 12 mCRC patient plasma samples (6 female, 6 male, average age 52 years, 52 draws) and 12 normal plasma (12 female) samples were purchased (ConversantBio, Huntsville, AL) and (ProMedix, Lafayette, LA) respectively. Five mCRC tissue samples were classified as KRAS mutation-positive by the vendor, but not tested as plasma. Additional samples (Fig.3) were provided by the Janku lab at MDACC. Samples were prepared using standard protocols (Qiagen CNA kit).
- ddPCR was performed on 1 – 8.75 μ l per sample per well using either the KRAS Screening Multiplex for ddPCR or validated PrimePCR ddPCR mutation assays for one of seven individual KRAS mutations (G12D, G12V, G13D, G12A, G12C, G12R, G12S, Bio-Rad Laboratories)
- Positive mutation references were from Horizon Diagnostics, and negative controls were wild type-only from Promega (female gDNA). Statistical significance was determined using 95% confidence intervals.

Figure 1: Multiplexed single-well detection of 7 actionable KRAS mutations




Figure 2: Cell-free plasma samples yield highly variable amounts of amplifiable DNA





Figure 3: KRAS Screening Multiplex assay on cell line DNA, FFPE, and cell-free plasma DNA



- A) 2D plots of KRAS Screening Multiplex assay on cell line DNA, FFPE, and cfDNA from MD Anderson Cancer Center patients.
- B) Concentration plot and Fractional Abundance plot (%mutant) of wildtype KRAS (green), mutant KRAS (blue), and % mutant (orange) from 6 MDACC cfDNA plasma samples, demonstrating range of %MT detection.

Figure 4: Multiplex detection of KRAS mutations in 24 patient samples from cell-free plasma DNA





Figure 5: ddPCR enables visualization of PCR inhibition from FFPE and cell-free plasma DNA samples



- 1D plots allow visualization and troubleshooting of PCR inhibition (A). For this sample (F3W), inhibitors are present (6 μ l loading), impacting positive fluorescence amplitudes. Loading less sample (2 μ l) allows better amplification. Regardless, the endpoint quantification is the same.

Conclusions

- Droplet digital PCR is an inexpensive method to quantify, absolutely, minimal amounts of FFPE and cfDNA, both for quantification and mutation detection.
- The amplifiable amount of cfDNA is significantly different between normal and cancer patient samples, and between KRAS mutant and KRAS wild type patient samples.
- We have demonstrated sensitive and precise detection down to 0.25% for multiple actionable KRAS mutations in cfDNA plasma samples from colorectal cancer patients.
- Droplet digital PCR provides a simple and robust workflow to screen a large volume of patient samples in a minimal amount of time.
- 1D & 2D plots of droplets enables rapid identification of PCR inhibition, either by poor assay design, sample inhibitors, poorly optimized conditions, or template degradation.

QX100 and QX200 Droplet Digital PCR Systems and/or their use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. Purchase of the product includes a limited, non-transferable right under such intellectual property for use of the product for internal research purposes only. No rights are granted for diagnostic uses. No rights are granted for use of the product for commercial applications of any kind, including but not limited to manufacturing, quality control, or commercial services, such as contract services or fee for services. Information concerning a license for such uses can be obtained from Bio-Rad Laboratories. It is the responsibility of the purchaser/end user to acquire any additional intellectual property rights that may be required.

